

UNITED STATES PATENT APPLICATION

FIDUCIAL MARKER DEVICES, TOOLS, AND METHODS

INVENTORS

Matthew S. Solar
of Indialantic, Florida, U.S.A.

Thomas L. Bridges
of Melbourne Beach, Florida, U.S.A.

David M. Lee
of Melbourne Beach, Florida, U.S.A.

Mark Freas
of Melbourne, Florida, U.S.A.

Attorney: Suneel Arora
Reg. No. 42,267
Schwegman, Lundberg, Woessner, & Kluth, P.A.
1600 TCF Tower
121 South Eighth Street
Minneapolis, Minnesota 55402
ATTORNEY DOCKET NO. 00723.077US1

resonance imaging (MRI), computed tomography (CT), positron emission tomography (PET), and single photon emission computed tomography (SPECT).

For example, in one technique, multiple fiducial markers are screwed into the patient's skull to define landmarks recognizable by an imaging system.

5 The imaging system is used to obtain one or more preoperative images of the patient's brain. Recognizable images of the fiducial markers appear on such preoperative images. Such a bone-anchored fiducial marker typically includes an externally threaded bone-screw portion, which is driven into the skull. A threaded shaft rises up and out of the skull from the bone-screw. The threaded
10 shaft typically receives a screwed-on imagable sphere that is visible on an MRI or CT image. The multiple fiducial markers on the patient's skull define landmarks on preoperative images that are useful to the physician for planning entry coordinates on the patient's skull and for planning a trajectory to a target location in the brain. An image-guided surgical workstation uses these
15 preoperative images and the planning data to guide the neurosurgeon while actually performing the subsequent surgical procedure.

After the preoperative planning phase, the patient is brought into the operating room so that the planned surgical procedure can be performed. On the operating table, the patient's skull is clamped in a head-frame or otherwise
20 immobilized. In order to use the preoperative images provided by the image-guided workstation to guide the surgeon during the surgical procedure, the patient's skull must first be "registered" to the preoperative images. The registration creates an association between (1) the actual physical location of the fiducial markers on the patient's skull in the operating room and (2) the locations
25 of the images of the fiducial markers visible on the preoperatively-obtained images. This allows mapping between the actual space in which the patient is located to the space defined by the preoperative images.

According to one registration technique, a "wand" is used to perform this patient registration. The wand typically includes multiple light-emitting diode
30 (LED) locators or reflective locators, which are visible to an infrared camera or other detector of an optical positioning system in the operating room. The camera and optical positioning system are operatively connected to the image-guided workstation. The locators define the position of the wand in the operating room, including the position of a sharp tip portion of the wand, which

is in a known physical relationship to the locators. To register the patient, the imagable spheres are unscrewed from the fiducial marker shafts, and replaced by respective "divots" that are sized and shaped to receive the wand tip. These divots are screwed or otherwise engaged onto the respective fiducial marker shafts, such that when the wand tip is received into the maximum depression point of the divot, the wand tip then corresponds to the same location as the center of the imagable sphere when the imagable sphere was screwed onto the fiducial marker shaft. A reference divot is typically also present in the operating room at a known location, such as attached to the operating table or the patient's skull-immobilizing head-frame. During the patient registration process, the surgeon touches the wand tip to the reference divot (to provide an absolute positional reference to the image-guided workstation), and then to each fiducial marker divot. This permits the image-guided workstation to correlate the actual physical location of the patient's skull to the preoperative images. The physician can then use the wand, in conjunction with the preoperative images provided by the image-guided workstation, to locate an appropriate entry point and trajectory to the target in the brain.

The present inventors have recognized that problems with the above registration procedure include patient discomfort caused by the presence of the fiducial markers, increased trauma to the patient resulting from using multiple fiducial markers screwed into different locations of the patient's skull, the difficulty of unscrewing the imaging spheres and replacing them with the registration divots, a limited field of view of the camera used in the operating room, and the difficulty of constructing a multi-modal fiducial marker that can be recognized by more than one imaging modality or positioning system. Moreover, the present inventors have recognized the desirability of streamlining the registration process to reduce its time and cost. For these and other reasons, which will become apparent upon reading the following detailed description and viewing the drawings that form a part thereof, the present inventors have recognized an unmet need for improved fiducial marker devices, tools, and methods.

BRIEF DESCRIPTION OF THE DRAWINGS

In the drawings, which are not necessarily drawn to scale, like numerals describe substantially similar components throughout the several views. Like numerals having different letter suffixes represent different instances of substantially similar components. The drawings illustrate generally, by way of example, but not by way of limitation, various embodiments discussed in the present document.

Figure 1A is a schematic diagram illustrating generally one example of an imagable fiducial marker that includes a built-in conical divot or other male or female receptacle, or the like.

Figure 1B is a schematic diagram illustrating generally one example of an imagable fiducial marker that omits the divot illustrated in Figure 1A, but which is both locatable by a remote positioning system and imagable by one or more imaging modalities.

Figure 2A is a schematic diagram illustrating generally an alternative example of a fiducial marker that includes a cylindrical imaging fiducial locator and a conical or other divot or other receptacle for receiving a positioning wand tip or the like.

Figure 2B is a schematic diagram illustrating generally one example of an imagable fiducial marker that omits the divot illustrated in Figure 2A, but which is both locatable by a remote positioning system and imagable by one or more imaging modalities.

Figure 3A is a schematic diagram illustrating generally one example of a positioning wand for use in conjunction with a remotely-located camera or other like device of an optical positioning system, such as can be coupled to an image-guided surgical workstation in an operating room.

Figure 3B is a schematic diagram, similar in certain respects to Figure 3A, illustrating generally one example of a positioning wand including energy reflective surfaces that are capable of being oriented or aimed toward a remote detector.

Figure 3C is a perspective view schematic diagram illustrating generally, by way of example, but not by way of limitation, certain generally “cylindrical” columnar structures having faceted lateral peripheral surfaces.

Figure 3D is a schematic diagram illustrating generally an example of a positioning wand with flat disk-shaped pieces of reflective tape are attached in a known configuration.

Figure 4 is a schematic diagram illustrating generally, by way of example, but not by way of limitation, an image guided surgical (IGS) computer workstation to which an optical positioning system is coupled.

Figure 5 is a schematic diagram illustrating generally a unitary divot assembly that includes multiple divots.

Figure 6A is a schematic diagram illustrating generally a divot assembly that includes a swiveling tilted head carrying a conical or other divot or the like.

Figure 6B is a schematic diagram illustrating generally a locator assembly that includes a swiveling tilted head including a surface that reflects electromagnetic energy.

Figure 7A is a schematic diagram illustrating generally a divot assembly that includes a swiveling and pivotable head carrying a conical or other divot.

Figure 7B is a schematic diagram illustrating generally a divot assembly that includes a swiveling and pivotable head including a surface that reflects electromagnetic energy.

Figure 8 is a schematic diagram illustrating conceptually a fiducial marker carrier that is attachable to (and also detachable from) a single location on the patient's skull, thereby reducing trauma to the patient.

Figure 9 is an exploded view schematic diagram illustrating generally one example of the carrier, including a frame, a post, and a base.

Figure 10 is a schematic diagram illustrating a portion of a fiducial marker carrier that includes at least one antirotational spike for engaging the surface of the skull.

FIG. 11 is a perspective view of an alternative example of a fiducial marker.

FIG. 12 is a top view of the fiducial marker illustrated in FIG. 11.

FIG. 13 is a perspective view of modified unitary fiducial marker.

FIG. 14 is a perspective view of an optional imagable plug.

FIG. 15 is a side view of the optional imagable plug of FIG. 14.

FIG. 16 is a perspective view of an optional fluid absorbing cover (or coating).

FIG. 17 is a side cross-sectional view of an alternative example of a fiducial marker.

FIG. 18 is a top view of an the fiducial marker of FIG. 17.

FIG. 19 illustrates a side view of a fiducial marker that includes a self-
5 drilling and self-tapping threaded distal tip portion.

FIG. 20 illustrates a side view of a fiducial marker that includes a threaded distal tip portion that need not be self-tapping and/or self-drilling.

FIG. 21 illustrates a side view of a fiducial marker that includes a barbed distal tip portion.

10 FIG. 22 illustrates a side view of a fiducial marker having a distal tip portion that includes tangs, or another laterally expandable retention element.

FIG. 23 illustrates a side view of a fiducial marker including a laterally expandable retention element and also having a self-tapping and/or self-drilling externally threaded distal tip portion.

15 FIG. 24 is a side cross-sectional view of a fiducial marker having a protective cap.

FIG. 25 is a side cross-sectional view of a protective cap with an adjustable-height skirt.

20 FIG. 26 is a side cross-sectional view of a protective cap disposed about a fiducial marker.

FIG. 27 is a perspective view of a protective collar that can be disposed about a fiducial marker that has been affixed to a subject's skull.

FIG. 28 is a perspective view of a protective collar and cap.

25 FIG. 29 is a perspective view of an alternate example of a protective collar that can be slipped over a fiducial marker.

FIG. 30 is a perspective view illustrating an example of a headband for protecting fiducial markers from mechanical impact.

FIG. 31 is a side view illustrating an example of a tubular protective guide collar.

30 FIG. 32 is a flow chart illustrating one example of using a guide collar.

FIG. 33 is a perspective view of an alternative guide collar.

FIG. 34 is a side sectional view of a guide base, a height extender, a fiducial marker, and a screwdriver shaft.

marker **100** that includes a built-in divot **102**. In this example, the divot **102** includes a female receptacle, such as the illustrated conical depression.

However, as used herein, a divot also refers to any other male or female receptacle, or the like. The divot **102** is capable of receiving a correspondingly sized and shaped mating tip of a positioning wand or like instrument. Such a wand or instrument is useful for registering the actual physical location of the patient's skull to preoperative or other images of the subject's brain. Such images are typically stored in a memory of an image-guided surgical (IGS) computer workstation.

10 In the example illustrated in Figure **1A**, the fiducial marker **100** includes an imagable substantially spherical fiducial locator **104**. The fiducial **104** is locatable using one or more imaging system modalities. In this example, a shaft **106** extends orthogonally outward from a circumferential portion of the spherical fiducial **104**. The shaft **106** includes an externally threaded portion **108**. The
15 externally threaded portion **108** is sized and shaped for being received within a correspondingly sized and shaped mating internally threaded receptacle **110** of an externally-threaded self-tapping base **112**. In this example, the base **112** is capable of being mounted in a skull **114**, such as either flush to (or even recessed from) an outer surface **116** of the skull **114**. One example of a suitable base **112**
20 is described in commonly-assigned Mazzocchi et al. U.S. Patent Application Serial No. 10/206,884 entitled FIDUCIAL MARKER DEVICES, TOOLS, AND METHODS, which was filed on July 24, 2002, and which is incorporated herein by reference in its entirety, including its disclosure relating to a flush or recessed mounted base and other fiducial marker devices, tools and methods. However,
25 in alternative examples, the base **112** need not be configured for mounting flush to or recessed from the outer surface **116** of the skull **114**. In this example, the shaft **106** includes a pointed tip **115**. This permits the shaft **106** to more easily penetrate a sterile drape that, in certain circumstances, may be placed over the patient's skull **114**. Moreover, in this example, the receptacle **110** of the base
30 **112** is shaped to accommodate the pointed tip **115**. However, in an alternative example, the tip **115** need not be pointed.

In one example, the imaging spherical fiducial locator **104** houses a generally spherical (e.g., except for the conic cutaway of the divot **102**) sealed interior cavity **118**. In one example, the cavity **118** is filled with an imagable

fluid that is visible on one or more imaging modalities (e.g., MR, CT, etc.). In this example, the apex of the conic divot **102** is located at a spherical center of mass of the imaging spherical fiducial locator **104** (i.e., the apex is located where the center of mass would be if the imaging fiducial locator **104** were perfectly spherical, without any cutout divot). This allows the tip of a positioning wand (recognizable by a camera in an optical position locating system that is coupled to the image-guided surgical workstation) to be inserted into the divot **102**. This results in the wand tip being located at the spherical center of mass of the imaging spherical fiducial locator **104**. This is useful for assisting in registering the physical location of the patient to the preoperative images stored in the image-guided surgical workstation.

Unlike fiducial marker assemblies that require the user to attach an imaging fiducial while obtaining the preoperative images of the patient's brain, and to then replace that imaging fiducial with a separate divot during patient registration in the operating room, the fiducial marker **100** illustrated in Figure **1A** does not require any such exchange of the imaging fiducial for a separate divot. Instead, the divot is integrated into the imaging fiducial itself, as illustrated in Figure **1A**. This reduces the complexity of the image-guided surgical procedure and, therefore, reduces its cost. It also reduces the complexity of manufacturing, which, in turn, reduces manufacturing costs.

In one example (but not by way of limitation), the base **112** is constructed of stainless steel. The shaft **106** and the imaging spherical fiducial locator **104** are constructed of molded plastic polymer. In this example, the imaging spherical fiducial locator **104** includes an open cavity **118** for receiving the imaging fluid, and for then receiving an insertable plastic conical divot **102** that adhesively or otherwise seals the cavity **118** to retain the imaging fluid therein. The imaging fluid in the cavity **118** is visible and provides good contrast on images produced by at least one imaging modality. In one example, the imaging fluid is multimodal (i.e., locatable by more than one imaging modality), such as by using a mixture of different imaging fluids that are locatable on different imaging modalities. In an alternative example, the plastic forming the imaging spherical fiducial locator **104** includes a substance that is viewable on a first imaging modality, while the imaging fluid within the cavity **118** is viewable on a different second imaging modality.

In one such illustrative example, the plastic imaging fiducial locator **104** is doped with a substance having a high atomic number (Z), such as barium, titanium, iodine, silver, gold, platinum, iodine, stainless steel, titanium dioxide, etc. that provide good contrast on a CT or other radiographic imaging system. In this illustrative example, the fluid within the cavity **118** includes gadopentatate dimeglumine, gadoteridol, ferric chloride, copper sulfate, or any other suitable MRI contrast agent, such as described in chapter 14 of Magnetic Resonance Imaging, 2nd ed., edited by Stark and Bradley, 1992, which is incorporated herein by reference.

In an alternative multimodal example, the cavity **118** is omitted. Instead, the spherical fiducial locator **104** is constructed of a substantially solid plastic or other material that is hygroscopic, that is, capable of receiving and retaining a fluid, such as an imaging fluid that is viewable on an imaging system (e.g., an MRI imaging system or the like). In a further example, the plastic forming the spherical fiducial locator **104** is doped or otherwise includes a substance that is viewable on a different imaging system, such as, for example, a CT or other radiographic imaging system. Illustrative examples of solid plastics that can be made hygroscopic include, among other things, nylon and polyurethane. Using a hygroscopic material avoids the complexity and cost associated with manufacturing a sealed cavity **118** for retaining an imaging fluid. Moreover, by adapting the solid hygroscopic plastic for imaging using a first modality, and by using the imaging fluid for imaging using a second modality, each of the solid and the fluid can be separately tailored toward providing better contrast for its particular imaging modality.

In another alternative example in which the cavity **118** is omitted, the fiducial locator **104** includes a rigid solid (e.g., substantially spherical, but for the conic divot) interior. This solid material is doped with a substance that provides good contrast using a first imaging modality (e.g., CT). A hygroscopic outer coating is formed thereupon. The coating permits soaking up a fluid that provides a good contrast using a second imaging modality (e.g., MRI).

In a further example of the fiducial marker **100** illustrated in Figure **1A**, the outer surface of the imaging spherical fiducial locator **104** is reflective of light or other electromagnetic energy. Consequently, it is also locatable by the operating room camera in an optical positioning system that is coupled to the

image-guided workstation (e.g., during patient registration). In one such example, the outer surface of the imaging spherical fiducial locator **104** includes light-reflective micro-spheres (e.g., embedded in an adhesive covering the imaging spherical fiducial **104**). In another such example, the outer surface of the imaging spherical fiducial **104** is covered with an adhesive-backed light-reflective tape, such as SCOTCHLITE® 9810 Reflective Material Multipurpose Tape sold by Minnesota Mining and Manufacturing Co. (“3M®”), of Saint Paul, Minnesota.

Figure **2A** is a schematic diagram illustrating generally, by way of example, but not by way of limitation, an alternative example of a fiducial marker **200** that includes a generally cylindrical imaging fiducial locator **202** and a conical or other divot **102**. In one example, the generally cylindrical imaging fiducial locator **202** includes a sealed cavity **204** for receiving and retaining an imagable fluid, as discussed above. In another example, the sealed cavity **204** is omitted, as discussed above. In one such example, the generally cylindrical imaging fiducial locator **202** is instead constructed of a substantially solid hygroscopic plastic that carries an imagable fluid (as discussed above), such as for providing multimodal contrast across different imaging modalities. In a further example, the generally cylindrical outer surface of the imaging fiducial locator **202** is reflective, as discussed above, such that the imaging fiducial locator **202** is also visible to a camera of an optical position locating system that is coupled to an image-guided surgical workstation (e.g., during patient registration and/or a subsequent image-guided surgical procedure). In one such example, the imaging fiducial locator **202** is covered with adhesive-backed reflective tape taken from a rectangular strip of such tape that is wound into a roll. In this example, the generally cylindrical shape of the outer surface of the imaging fiducial locator **202** is much easier to wrap using a wound rectangular strip of the adhesive reflective tape than a spherical surface, such as is illustrated in Figure **1A**, and therefore costs less to manufacture. In this document, the term “generally cylindrical” is not limited to a perfectly cylindrical surface, but instead is understood to include any faceted or other column or like structure (e.g., an octagonal cylinder a hexagonal cylinder, etc.) that includes a lateral peripheral surface that easily accommodates receiving a wound rectangular or similar strip of tape (as opposed to a spherical, elliptical, or conical surface, to

which is more difficult to evenly apply a wound rectangular strip of tape taken from a roll). Examples of such generally “cylindrical” columnar structures having faceted lateral peripheral surfaces are illustrated in Figure 3C.

In an alternate example to the illustrations of Figures 1A and 2A, the
5 divot **102** is omitted from the fiducial marker **100** or **200**. However, the
resulting fiducial marker is still configured to be locatable by a remote
positioning system as well as imagable using one or more imaging modalities.
In one such example, the outer surface **104** or **202** is still configured to be light
reflective, such as discussed above. In one such example, the fiducial markers
10 **100** and **200** still advantageously are locatable using one or more imaging
modalities (e.g., MR, CT, or other imaging system providing 3D or other
internal images within a subject) as well as also being locatable external to the
subject, such as by using a remote camera or like component of an optical or
other positioning system, e.g., that is coupled to an image-guided workstation.
15 In one example, this permits automatic registration of the actual location of the
subject in the operating room (e.g., using the cameras to locate the light
reflective fiducial markers **100** or **200**) to preoperative images of the patient on
which the same imagable fiducial markers **100** and **200** appear. This eliminates
any need to register the patient by inserting an optically-locatable positioning
20 wand tip into a divot of each fiducial marker (and also eliminates any need for a
reference divot or other absolute position reference), because the fiducial
markers themselves are optically locatable and registerable to known locations
on the preoperative images. Therefore, in this example, the divots **102** are not
needed and can be omitted, as illustrated by the divotless spherical imagable
25 reflective fiducial marker **120** in Figure 1B and the divotless cylindrical
imagable reflective fiducial marker **206** in Figure 2B. Although Figure 2B
illustrates an example including a cavity **204** for carrying a liquid contrast agent,
in an alternative example, the cavity **204** is omitted, and the fiducial marker **206**
includes a solid structure that is doped or otherwise configured (e.g.,
30 hygroscopic) for providing good imaging contrast using one (e.g., CT) or more
imaging modalities.

In yet another example, the fiducial markers **100** and **200** respectively
illustrated in Figures 1A and 2A include the illustrated divots **102** and are
locatable by a remote positioning system (such as by including light-reflective

outer surfaces and/or embedded coils that perform magnetic field sensing in a magnetic field based positioning system). However, in this example, the fiducial markers **100** and **200** need not be configured for providing contrast on the one or more imaging modalities. In such an example, the preoperative images are taken
5 with imagable fiducial markers placed within respective bases **112**. Such imagable fiducial markers are then replaced (within their respective bases **112**) by nonimagable fiducial markers that are locatable by a remote positioning system, such as by including both a divot and a light-reflective surface. The light reflective surface permits automatic location by the remote positioning
10 system. However, if the reflective surface is dirty or otherwise unrecognizable by the remote positioning system, a wand or other locating instrument can be placed within the divot to perform the remote locating of the fiducial marker.

Moreover, although Figures **1A** and **2A** illustrate examples in which a shaft **106** is received within a base **112** that is mounted flush to (or recessed
15 from) the outer surface **116** of the skull **114**, this is not required. In one alternate example, the shaft **106** is manufactured as a stainless steel or other suitable material that is capable of acting as a self-tapping bone screw. In such an example, the threaded portion **108** of the shaft **106** is threaded directly into the skull **114** without using any base **112**. In another alternate example, the base
20 **112** includes a shaft or flange portion that rises above the outer surface **116** of the skull **114**. In certain examples, the fiducial markers **100** and **200** may use a threaded or other shaft **106** for coupling to the base **112**, or alternatively may use a snap-fit clip or a like attachment device for coupling to the base **112**.

Figure **3A** is a schematic diagram illustrating generally, by way of
25 example, but not by way of limitation, one example of a positioning wand **300**, such as for use with a remotely-located camera or other like device of an optical positioning system configured for being coupled to an image-guided surgical workstation in an operating room. In this example, the wand **300** includes a tip
30 **302** that is sized and shaped to permit being received in a divot **102** of a skull-mounted fiducial marker (such as fiducial markers **100** and **200**). The wand **300** includes a plurality of cylindrically-shaped fiducial locators **304** that are locatable by the camera or other like device of the optical positioning system. The fiducial locators **304** (which typically need not include divots) on the wand **300** are positioned in a known spatial relationship to each other and to the tip

302 of the wand 300. By recognizing the locations of the fiducial locators 304, the optical positioning system is capable of computing the location of the wand tip 302, which is in a known spatial relationship with the configuration of fiducial locators 304. This permits the wand 300 to be used in conjunction with the optical positioning system to register the patient and to further plan and/or perform the surgical procedure using the image-guided surgical workstation. The fiducial locators 304 are covered with adhesive-backed reflective tape, as discussed above. The cylindrical (or faceted cylindrical) shape of the fiducial locators 304 permits easier wrapping by the reflective tape than the spherical fiducials, as discussed above. This reduces the cost of manufacturing the fiducial locators 304 and, in turn, reduces the cost of manufacturing the positioning wand 300.

Figure 3B is a schematic diagram, similar in certain respects to Figure 3A, but illustrating a wand 306 that includes locators 308A-C having swiveling or fixed cylindrical locators 308A-C having respective slanted (e.g., flat, parabolic, or other) top surfaces 310A-C (e.g., non-orthogonal with respect to a longitudinal center axis 311 of the locator 308) that reflect light or other electromagnetic energy for being located by a remote detector. In an example in which the locators 308A-C swivel, each such locator 308 includes a shaft inserted into a hole or other receptacle in the wand 306. This permits the locator 308 to rotate with respect to its mounting location on the wand 306. Either the wand 306 itself or the individual locators 308A-C are oriented by the user to aim the reflective surfaces 310A-C toward a camera or other detector of an optical positioning system. In one further example, the circumferential surfaces of the cylindrical locators 308A-C are also light-reflective, however, this is not required. In one such cost-effective example, the reflective tape disks are adhered to the flat slanted top surfaces 310A-C and the circumferential lateral surfaces of the cylindrical locators 308A-C are not reflective.

Figure 3C is a perspective view schematic diagram illustrating generally, by way of example, but not by way of limitation, certain generally “cylindrical” columnar structures 312, 314, and 316 having faceted lateral peripheral surfaces. Such surfaces are conducive to receiving a rectangular or like strip of adhesive reflective tape. Such structures, therefore, are particularly well-suited for implementing locators that are remotely locatable by an optical positioning

system. Such remotely detectable locators are suitable for use in the fiducial markers illustrated in Figures 2A and 2B, as well as for use in the remotely detectable locators of the positioning wands illustrated in Figures 3A and 3B. Such remotely detectable locators are also useful for being affixed in a known relationship to the patient, such as to the operating table or to a skull-immobilizing headframe. This provides a remotely detectable absolute positional reference to an optical positioning system. Such remotely detectable locators are also useful for being affixed to a biopsy needle, shunt catheter, or other instrument being introduced through a trajectory guide device or otherwise used in an image-guided surgical procedure.

Figure 3D is a schematic diagram illustrating generally, by way of example, but not by way of limitation, an alternative example of a positioning wand 318. In this example, which flat disk-shaped pieces of reflective tape are attached to the wand 318 in a known configuration, such as at the distal ends of radial arms extending therefrom.

Example 2

Figure 4 is a schematic diagram illustrating generally, by way of example, but not by way of limitation, an image guided surgical (IGS) computer workstation 400, which is capable of displaying previously acquired and loaded preoperative images of a patient's skull. On these preoperative images appear viewable images of imagable fiducial markers that were screwed into the patient's skull before the preoperative imaging (e.g., using MRI, CT, etc.). In the example illustrated in Figure 4, the imagable fiducial locators have been unscrewed from respective bases 402 screwed into the patient's skull. The imagable fiducial locators have been replaced by patient registration divot assemblies 404 that have been screwed into (or otherwise coupled to) respective bases 402 in the patient's skull 114. In this example, the registration divot assemblies 404 are configured to receive a shaft tip 406 of a positioning wand 408 that is locatable by one or more remote cameras 410A-B (or other sensing devices) of an optical position detection system 412 connected to the IGS workstation 400. In one example, the positioning wand 408 includes spherical reflective fiducial locators 414. The fiducial locators 414 are arranged in a known spatial relationship to each other (however, it may alternatively use other reflective locators such as discussed elsewhere in this document). The optical

positioning system 412 includes an infrared light (or other energy source) 416 that provides light that is reflected from the reflective fiducial locators 414. This permits the reflective fiducial locators 414 on the positioning wand 408 to be located and recognized by the cameras 410A-B. In some circumstances, however, the field of view (or “sweet spot” of the field of view) provided by cameras 410A-B is limited. This sometimes makes it difficult for the optical positioning system 412 to recognize the positioning wand 408. Moreover, the recessed receptacle in the divot assembly 404 typically limits the range within which the probe 408 can be manipulated (e.g., to bring it within the field of view) while retaining the wand tip 406 within the recessed receptacle.

Figure 5 is a schematic diagram illustrating generally, by way of example, but not by way of limitation, a unitary divot assembly 500 that includes multiple divots 502. In this example, the unitary divot assembly 500 is configured such that it can be threaded into or otherwise coupled to a base 504 that is secured to the patient’s anatomy (wherein the base 504 is also configured for alternatively receiving an imagable fiducial locator, e.g., during preoperative imaging). Figure 5 illustrates multiple conical receptacle divots 502 having commonly located apexes. These commonly located apexes are designed to coincide with the center of the image produced by the imagable fiducial locator for which the divot assembly 500 has been substituted during patient registration. In the illustrated example, the divots include a top conical divot 502A and four side conical divots 502B-F. The four side conical divots 502B-F are distributed around the cylindrical lateral peripheral circumference of the upper portion of the divot assembly 500. The wand tip 406 may be inserted into any one of the divots 502. This permits a greater range of motion of the positioning wand 408. As a result, it is easier to bring the reflective fiducials 414 on the positioning wand 408 into the field of view of the cameras 410A-B of the optical positioning system 412.

Figure 6A is a schematic diagram illustrating generally, by way of example, but not by way of limitation, a divot assembly 600 that includes a swiveling tilted head 602 carrying a conical or other divot 604 or the like. In this example, the head 602 is tilted with respect to a cylindrical coupling 606 extending outwardly therefrom. The coupling 606 includes a hollow interior or other (female or male) connector that snap-fits onto and rotatably rides upon a

5 mating (male or female) connector 608 that is located at a proximal end of a shaft 610 portion of the divot assembly 600. The swiveling apex 612 of the divot 604 is designed to coincide with the center of mass of the imagable fiducial locator for which the divot assembly 600 has been substituted during patient registration. The swiveling tilted head 602 permits a wide range of motion of the positioning wand 408 when the wand tip 406 is inserted into the divot 604. As a result of such rotational articulation, it is easier to bring the reflective fiducial locators 414 on the positioning wand 408 into the limited field of view of the cameras 410A-B of the optical positioning system 412.

10 Figure 7A is a schematic diagram illustrating generally, by way of example, but not by way of limitation, a divot assembly 700 that includes a swiveling and pivotable head 702 carrying a conical or other divot 704. In this example, the head 702 is carried by a shackle-like U-shaped bracket 704 that rotatably rides upon a snap-fit or other capturing post 706 that extends upward
15 from a shaft portion 708 of the divot assembly 700. This allows swiveling of the bracket 704 (and the head 702 carried by the bracket 702) with respect to the shaft 708. In this example, the head 702 is suspended between upward-projecting risers of the bracket 704 by axels 710A-B extending outward from opposing sides of the head 702 and received within corresponding receptacles in
20 the risers of the bracket 704. This permits pivoting/tilting articulation of the head 702 with respect to the swiveling bracket 704. Therefore, this example provides a swiveling and adjustably tiltable divot 704 that is designed such that its apex 712 coincides with the center of mass of the imagable fiducial locator for which the divot assembly 700 has been substituted during patient
25 registration. Among other things, the swiveling tiltable head 702 advantageously permits a greater range of motion of the positioning wand 408 when the wand tip 406 is inserted into the divot 704. As a result, it is easier to bring the reflective fiducials 414 on the positioning wand 408 into the limited field of view of the cameras 410A-B of the optical positioning system 412.

30 Figures 6B and 7B are schematic diagrams that are similar in certain respects to Figures 6A and 7A. However, the locator assemblies 614 and 714 illustrated by respective Figures 6B and 7B omit the respective divots 604 and 704. Instead, the locator assemblies 614 and 714 provide aimable electromagnetic energy (e.g., light) reflective surfaces 616 and 716, respectively.

The reflective surfaces **616** and **716** are aimed at the camera of an optical positioning system **412** to allow automatic detection of the locator assemblies **614** and **714** without requiring the use of a positioning wand **408**.

5 The reflective surfaces **616** and **716** are configured so that, when aimed properly, they produce a reflected image that can be correlated to a previously acquired patient image on which an image of an imagable fiducial marker appears. In one such example, reflective surface **616** corresponds to the center of mass of a similarly sized spherical locator on an imagable fiducial marker assembly for which locator assembly **614** is substituted during patient
10 registration. In another such example, reflective surface **716** includes a circular disk-shaped piece of reflective tape affixed to a surface **718** such that this reflective disk pivots about the axis provided by axels **710A-B**. In this manner, the reflected disk shape corresponds to the center of mass of a similarly sized spherical locator on an imagable fiducial marker assembly for which locator
15 assembly **714** is substituted during patient registration.

Example 3

As discussed above, screwing multiple fiducial markers into different locations in the patient's skull **114** results in trauma and/or risk of infection at each one of such multiple different locations. Figure **8** is a schematic diagram
20 illustrating conceptually, by way of example, but not by way of limitation, a fiducial marker carrier **800** that is attachable to (and also detachable from) a single location on the patient's skull **114**, thereby reducing trauma and risk of infection to the patient. In this example, the fiducial marker carrier **800** is configured for carrying multiple different imagable fiducial locators **802** such
25 that they are positioned at different locations about the patient's skull **114**. As discussed below, the carrier **800** uses a keyed mounting arrangement, such that the carrier **800** can be attached to the patient's skull **114**, then detached from the patient's skull **114**, and later reattached to the patient's skull **114** in the same orientation in which it was initially attached to the patient's skull **114**.

30 In the example illustrated in Figure **8**, the carrier **800** includes a keyed frame **804** that is attached to a keyed post **806** for mounting. The keyed post **806** is, in turn, attached to a single flush-mounted or recessed-mounted or other keyed base **808**, which was previously screwed into the patient's skull **114**. This keyed arrangement of the frame **804**, the post **806**, and the base **808** permits

attachment, detachment, and reattachment in the same orientation as the original attachment, as discussed above. In an alternative example, the post 806 is integrally formed as part of the frame 804, rather than being keyed for attachment thereto.

5 In one example, such illustrated in Figure 8, the imagable locators 802 are placed about the subject's head such that they surround the patient's skull. Although such a surrounding arrangement is not required, it is believed to improve the accuracy of using the images of the locators 802 (e.g., in conjunction with the IGS workstation) for planning and/or performing an image-guided surgical procedure, as compared to an arrangements in which locators are
10 disposed more closely together (e.g., on the same side of the subject's head).

 Figure 9 is an exploded view schematic diagram illustrating generally, by way of example, but not by way of limitation, one example of the carrier 800, including the frame 804, the post 806, and the base 808. In this example, the
15 base 808 includes self-tapping external threads 902, and is capable of being mounted flush with (or even recessed within) the patient's skull 114. The base 808 includes an internally-threaded receptacle 904 that is sized and otherwise configured such that it is capable of receiving a screw. The base 808 also includes a female or male keying feature for receiving a mating keying feature of
20 the post 806 to fixedly define the orientation of the post 806 with respect to the base 808. In one example, the keying feature includes a key slot 906 extending radially outward from the receptacle 904 along a proximal surface of the base 808.

 The post 806 includes a proximal end 908 and a distal end 910. The post
25 806 includes a center lumen 912 in which an attachment screw 914 is received and seated. The screw 914 attaches the post 806 to the base 808. The distal end 910 of the post 806 includes a male or female keying feature (such as a key protrusion 916 extending radially outward from the center lumen 912 along the distal end 910 of the post 806) that mates with the keying feature (e.g., key slot
30 906) of the base 808. Such mating during the attachment fixedly defines the orientation of the post 806 with respect to the base 808.

 In this example, the center lumen 912 includes a keyed seating receptacle 918 (or an analogous male keyed feature) for receiving a mating keyed feature of the frame 804. In the illustrated example of Figure 9, the keyed seating

receptacle 918 includes an increased diameter of the center lumen 912 (with respect to more distal portions of the center lumen 912) to provide the seating, and a radially-outwardly extending slot 920 to provide the keying.

In the example illustrated in Figure 9, the frame 804 includes legs 922A-D (or a fewer or greater number of legs 922), such as extending radially outwardly from a hub 924 and downwardly toward the middle portion of the patient's skull. Each of the legs 922 includes, such as at its respective distal end, a threaded receptacle 924A-D (or a snap-fitting or any other coupling) for receiving at least one of an imagable fiducial marker assembly 926, a divot assembly 928, a locator assembly 930 (e.g., reflector, LED, microcoil, etc.) that is remotely detectable by a positioning system in an operating room, or a combination 932 of two or more of the above. In an alternative embodiment (for example where a combination 932 includes an imagable locator and at least one of an operating room position locator and a divot), instances of such a combination 932 may be permanently affixed to corresponding locations on the legs 922 of the frame 804.

In the example illustrated in Figure 9, the hub 924 portion of the frame 804 also includes a downwardly protruding key 934 (or analogous female receptacle) that mates to the keyed seating receptacle 918, of the post 806, into which the key 934 is received. This fixedly defines the orientation of the frame 904 with respect to the post 806. A screw 936 is inserted through the hub 924, the key 934, and into an engaging interior threaded portion of the center lumen 912. This securely attaches the frame 904 to the post 806 in the fixedly defined orientation. The example illustrated in Figure 9 also includes at least one optional instrument mount 938. In one example, a reference divot (e.g., providing a position reference) is attached to the instrument mount 938.

Although Figures 8 and 9 illustrate examples in which a fiducial marker carrier 800 is mounted using a single base 808, in other examples, the carrier may be mounted using two or more bases 808 at the same location on the patient's skull (that is, at adjacent locations within the same scalp incision, or like limited trauma/infection risk zone; the incision need only be large enough to accommodate the two or more bases 808). Using two or more side-by-side bases 808 to attach the post 806 avoids potential rotational misalignment of a single base 808 coming slightly unscrewed from its original position.

Alternatively, if a single base **808** is used, such rotational misalignment can be avoided by including one or more antirotation spikes **1000** on the bottom of the distal end **910** of the post **806**, such as illustrated generally in Figure 10. In the example illustrated in Figure 10, the distal end **910** of the post **806** is keyed both to the base **808** and, using the antirotation spike(s) **1000**, to indentation(s) made in the surface **116** of the skull **114**. However, in an alternative example, the post **806** and the base **808** need not be keyed to each other. Instead, in such an example, the post **806** is keyed only to indentation(s) made by the antirotation spike(s) **1000** in the surface **116** of the skull **114**.

Example 4

FIG. 11 is a perspective view of an alternative example of a fiducial marker **1100**. FIG. 12 is a top view of the fiducial marker **1100** illustrated in FIG. 11. In the example of FIGS. 11 - 12, the unitary fiducial marker **1100** includes a substantially spherical head **1102**. A unitary fiducial marker includes both a single piece as well as multiple pieces that are assembled into a single assembly that, in use, is not disassembled or otherwise decomposed into more than one separate component. In this example, a divot **1104** is cut out from a proximal portion of the head **1102**. The divot **1104** is shaped to receive a corresponding mating shaped portion of a remote positioning locator. In one illustrative example, the divot **1104** is conical (as illustrated in FIG. 11), such as to receive a mating conical tip **302** of the positioning wand **300** illustrated in FIG. 3, or a similar probe tip. An apex of the inverted conical divot **1104** corresponds to a centroid of the substantially spherical head **1102**. In this example, a bone screw shaft **1106** extends outward from an opposite (e.g., distal) portion of the head **1102**. (Alternatively, if a sterile drape or the like is to be used between the tip **302** of the wand **300** and the divot **1104** of the fiducial marker **1100**, then, in one example, the location of the apex of the divot **1104** may be adjusted to offset the thickness of the sterile drape such that the tip **302** of the wand **300** is located at the centroid of the head **1102** even when the drape is interposed between the tip **302** and the divot **1104**).

In this example, the conical divot **1104** of the head **1102** includes slots **1108** extending therefrom. The slots **1108** accommodate a driving tip of a screwdriver (e.g., Phillips and/or flathead, etc.). In this manner, the slots **1108** permit the fiducial marker **1100** to be screwed into a skull, bone, or other

structure. Alternatively, the divot **1104** includes any other known rotational engagement structure for permitting rotation of the fiducial marker **1100** for threading it into bone, as discussed below.

In one example, the shaft **1106** includes one or more self-tapping or other external bone screw threads **1110**, which are sized and shaped for being threaded into bone, such as a patient's skull. In one example, a distal tip of the shaft **1106** includes at least one cutout, such as a quarter cylindrical cutout **1112**. In this example, the vertically-oriented flute-like cutout **1112** portion of the shaft **1106** assists in cutting bone as the shaft **1106** is being turned for threading into the bone. The self-drilling cutout **1112** and self-tapping nature of the threads **1110** are not essential. These features are not needed, for example, where a pre-drilled hole is available and used for receiving the shaft **1106**.

In one example, the unitary fiducial marker **1100** is made from substantially pure or alloyed titanium, substantially pure or alloyed stainless steel, and/or a ceramic. In one example, the resulting substantially spherical head **1102** is radiolucent and/or radiographically imagable and viewable using computed tomography (CT).

In the example of FIG. 11, the unitary fiducial marker **1100** includes an imagable locator head **1102** that is spherical (or otherwise shaped) for obtaining accurate location information (e.g., of its center). The head **1102** also includes a receptacle (such as the divot **1104**) that is shaped for receiving a mating portion (e.g., tip **302**) of a positioning instrument (e.g., wand **300**) during patient registration. Therefore, the unitary fiducial marker **1100** (with integrated imaging and registration divot) in the example of FIG. 11 avoids having to replace an imagable portion of a two-piece fiducial marker (used during preoperative imaging) with a separate registration divot (used during patient registration in the operating room). This simplifies an image-guided surgical procedure using the unitary fiducial marker **1100** having both the imagable head **1102** and the integrated divot **1104**. Such simplification should help lower the cost of the image-guided surgical procedure.

FIG. 13 is a perspective view of modified unitary fiducial marker **1100**. In this example, the shaft **1106** includes a threaded distal portion **1300** and an unthreaded proximal portion **1302**. The unthreaded proximal portion **1302** distances the head **1102** from the surface into which the threaded distal portion

1300 is screwed. In this example, the unthreaded proximal portion 1302 of the shaft 1106 is of a larger cylindrical diameter than the tapered threaded distal threaded portion 1300 of the shaft 1106. This forms a circular shoulder or seat 1304 at the base of the unthreaded proximal portion 1302 where it meets the threaded distal portion 1300. When the seat 1304 is of a larger diameter than the major diameter of the threads 1110, the seat 1304 provides a shoulder acting as a depth stop that inhibits the fiducial marker 1100 from being further advanced into the bone, such as by an accidental impact to the head 1102 of the fiducial marker 1100 that produces a mechanical shock.

10 In one example, several fiducial markers 1100 are packaged and sold together as a kit. In one such example, such a kit includes two or more different fiducial markers 1100 having different lengths of the unthreaded proximal portion 1302 of their respective shafts 1106. This accommodates patients having different skin or scalp thicknesses. For example, it may be desirable to keep the head 1102 portion of the fiducial marker 1100 above the patient's skin or scalp, while remaining as close to the skull as possible. If this is desired, it can be accomplished by selecting from the kit a particular fiducial marker 1100 having an appropriate shaft 1106 length to accommodate the skin or scalp thickness of the patient.

15 In this example, the seat 1304 includes a circular groove, channel, or kerf 1306. In this example, the kerf 1306 extends along the seat 1304 circumferentially around the threaded distal portion 1300. The kerf 1306 accommodates therein loose bone fragments that are channeled upward by the threads 1110 when the fiducial marker 1100 is being screwed into the skull.

20 Such groove, channel, or kerf 1306 for accommodating channeled bone fragments could similarly be incorporated into a distal side of the head 1102 in the examples of FIGS. 10 – 11, in which the threaded portion of the shaft 1106 extends directly from the head 1102.

Example 5

30 FIG. 14 is a perspective view of an optional imagable plug 1400. FIG. 15 is a side view of the optional imagable plug 1400, which can be made from the same material as the head 1102, if desired. The imagable plug 1400 is sized and shaped to be inserted into the divot 1104 during imaging such that the head 1102 presents a uniformly shaped imagable sphere to the imaging modality.

This assists in easier location of the centroid of the spherical combination of the head 1102 and the plug 1400, but is not believed to be required. In this example, the imagable plug 1400 is then removed during registration, thereby permitting access to the divot 1104. In one example, the plug 1400 includes fins 1402 that
5 are sized and shaped for engaging the corresponding slots 1108. In an alternative example, however, the fins 1402 are omitted.

In an alternative example, the imagable plug 1400 is made from a material having a slightly or substantially different imaging contrast property from the material comprising the rest of the head 1102. In this manner, an image
10 of the fiducial marker can be obtained in which the divot 1104 appears with a different imaging contrast than the rest of the head 1102. This shows the user where the divot 1104 is located within the image.

Example 6

FIG. 16 is a perspective view of an optional hydrophilic or hygroscopic
15 foam or other magnetic resonance (MR) imagable cover 1600 for slipping over the substantially spherical head 1102. In this example, the fluid/gel-carrying, fluid/gel-absorbing, or other fluid/gel-incorporating cover includes a circular or similar opening 1602 permitting the shaft 1106 to extend therethrough. In one example, a sterile and biologically safe magnetic resonance (MR) imagable
20 fluid/gel is soaked into the cover 1600 either before or after it is slipped over the head 1102. This allows the head 1102 to be imaged by MR as well as CT. In an alternative example, such multi-modality of imaging is similarly implemented using a preformed MR-imagable or other coating upon the head 1102, thereby avoiding any need for slipping a separate cover 1600 over the head 1102. Such a
25 fluid/gel-carrying, fluid/gel-absorbing, fluid/gel-incorporating, or other MR-imagable or other coating could be formed on the external spherical portion of the head 1102, or could additionally be formed in the divot 1104 as well. Examples of suitable coatings capable of soaking up an MR-imagable fluid or gel include, by way of example, but not by way of limitation: foam, silicone,
30 etc. Examples of MR imagable fluids for soaking into the cover 1600 (or coating) include, by way of example, but not by way of limitation: sterile saline, sterile saline or another fluid or gel mixed with gadolinium or another MR-imaging enhancing substance, etc.

Example 7

FIGS. 17 – 18 are side cross-sectional and top views, respectively, of an alternative example of a fiducial marker 1700 that is similar in certain respects to the example of FIG. 13. In FIGS. 17 – 18, the fiducial marker 1700 includes a substantially spherical head 1702. The head 1702 includes a conical or other
5 divot 1704 at its proximal side, and a shaft 1706 extending outwardly from its distal side. In this example, the shaft 1706 includes a proximal portion 1708 and a threaded distal tip portion 1710. The proximal portion 1708 and the threaded distal tip portion 1710 are separated by a shoulder or other seat 1712, such as described above. In this example, the divot 1704 of the head 1702 includes
10 rotational engagement features, such as slots 1714, for receiving a Phillips and/or flathead screwdriver or other driver. Alternatively, an Allen-type receptacle, or any other rotational engagement feature could be used for receiving another driver.

In the example of FIGS. 17 – 18, the head 1702 is made of a different
15 material than the shaft 1706. In one example, the different materials are selected to provide different image contrasts on a particular imaging modality (e.g., an MR image, a CT image, or even both types of images). In one such example, the head 1702 is relatively more highly visible on the particular imaging modality, and the shaft 1706 is less highly visible on the particular imaging modality.

20 In one example, this is effected by using a titanium shaft 1706 that includes a proximally projecting post 1716. In one example, a proximal end of the post 1716 provides the slots 1714, as illustrated in FIG. 17. In another example, the slots are instead incorporated into the head 1702. In this example, the head 1702 is a plastic sphere-like object that is insert-molded or otherwise
25 formed about the post 1716. In one example, the external surface of the post 1716 is knurled or roughened to promote adhesion of the head 1702 to the post 1716, such as during the insert-molding process. In one example, the head 1702 is highly MR-visible, while the shaft 1706 is not so highly MR-visible, but instead is radiolucent. In addition to insert-molding, other techniques for
30 affixing the head 1702 to the shaft 1706 include, without limitation, gluing, casting, spin-welding, and ultrasonic welding. In yet another example, the post 1716 is threaded, and the head 1702 is threaded and glued onto the post 1716.

Example 8

FIGS. 19 – 23 illustrate various distal tip configurations and techniques of attaching fiducial markers to bone. FIG. 19 illustrates a side view of a fiducial marker **1900** that includes a self-drilling and self-tapping threaded distal tip portion **1902**. This example may additionally include a vertical flute-like cutout, as discussed above, for enhancing its self-drilling capability. The head **1903** of the fiducial marker **1902** includes a conical or other divot **1904** and associated slots **1906** or other rotational engagement features for driving the fiducial marker **1900** into bone.

FIG. 20 illustrates a side view of a fiducial marker **2000** that includes a threaded distal tip portion **2002** that need not be self-tapping and/or self-drilling, such as for use when a hole has been pre-drilled into bone for receiving the tip portion **2002**. In one such example, the distal tip portion **2002** is neither self-tapping, nor self-drilling. In another such example, the distal tip portion **2002** is self-tapping, but is not self-drilling.

FIG. 21 illustrates a side view of a fiducial marker **2100** that includes a barbed or other distal tip portion **2102** enabling the fiducial marker **2100** to be driven into bone like a nail or a staple—that is, without needing any rotation. In one example, barbs **2104** help retain the distal tip portion **2102** within the bone. In another example, a nail-like distal tip portion **2102** is used instead. The nail-like distal tip portion **2102** may include a faceted point. In another example, the nail-like distal tip portion **2102** includes anti-rotation features that do not substantially inhibit the distal tip portion **2102** from being driven into bone, but which inhibit rotation after the distal tip portion **2102** has been driven into bone. The fiducial marker **2100** may be removed by grasping and pulling the proximal head **1903**, such as with a staple-puller-like tool. Therefore, this example need not include the slots **1906** or other rotational engagement features because rotation is not needed for inserting or removing the fiducial marker **2100**.

FIG. 22 illustrates a side view of a fiducial marker **2200** having a distal tip portion **2202** that includes tangs **2204A-B**, or another laterally expandable retention element. In one example, the tangs **2204A-B** are pushed outward by an ascending and/or descending longitudinally extending internal rod **2206** that pushes upward or downward against tapered internal shoulders of each of the tangs **2204A-B**. This, in turn, pushes the tangs **2204A-B** laterally outward in opposite directions. The rod **2206** extends longitudinally through an interior

passage **2208** of a shaft **2210**. The shaft **2210** extends between the distal tip **2202** and a head **2212** portion of the fiducial marker **2200**. In one example, the rod **2206** terminates at a proximal externally threaded drive head **2214** that engages an internally threaded portion of the head **2212**. The drive head **2214** includes screwdriver slots or one or more other rotational engagement features for turning the drive head **2214**. In one example, turning the drive head **2214** in a clockwise direction moves the drive head **2214** closer to the distal tip **2202** of the fiducial marker **2200**. This pushes the rod **2206** downward, which, in turn, pushes the tangs **2204A-B** outward to grip bone surrounding a pre-drilled hole into which the distal tip **2202** has been inserted. In another example, turning the drive head **2214** in a counter-clockwise direction moves the drive head **2214** away from the distal tip **2202** of the fiducial marker **2200**. This pulls the rod **2206** upward, which, in turn, pushes the tangs **2204A-B** outward to grip bone surrounding a pre-drilled hole into which the distal tip **2202** has been inserted.

FIG. **23** illustrates a side view of a fiducial marker **2300**, similar to the fiducial marker **2200** of FIG. **22**, but having a self-tapping and/or self-drilling externally threaded distal tip portion **2302**, such as for being introduced into bone without using a pre-drilled hole. The head **2212** of the fiducial marker **2300** of FIG. **23** also includes slots **1906** or other rotational engagement features for rotationally driving the fiducial marker **2300** into bone, such as by using a screwdriver. Then, the tangs **2204A-B** are forced outward as described above with respect to the fiducial marker **2200** of FIG.

Example 9

After a fiducial marker has been introduced into a patient's skull or other bone, it may be desirable to protect the fiducial marker, such as against accidental shocks or impacts, "twiddling" by the patient, etc.

FIG. **24** is a side cross-sectional view of a fiducial marker **2400** having a substantially spherical head **2402** that includes an internally threaded proximal divot **2406**, and a shaft **2408** extending outward from a distal side of the head **2402** toward a distal tip **2410** that has been threaded into a portion of the subject's skull **2412**. In this example, a protective cap **2414** has been threaded into the divot **2406**. The protective cap **2414** includes a disk-like top portion **2416** and a cylindrical circumferential skirt **2418**.

FIG. 25 is a side cross-sectional view of a further example of the protective cap 2414 in which the skirt 2418 includes an adjustable height outer cylindrical circumferential skirt 2420. In this further example, threads on the internal portion of the skirt 2418 engage threads on the outer portion of the skirt 2420, providing height adjustability to accommodate different scalp thicknesses. In use, the fiducial marker 2400 is first affixed to the subject's skull, then the protective cap is threaded into the divot 2406, and then the outer skirt 2420 is lowered to the appropriate height for the particular patient's scalp thickness.

FIG. 26 is a side cross-sectional view of another example of a protective cap 2600, which is disposed about a fiducial marker 2602 that has been affixed to a subject's skull. In this example, the cap 2600 includes a proximal disk portion 2604, a cylindrical circumferential portion 2606, and a distal base ring flange portion 2608. The distal base ring flange portion 2608 includes a self-adhesive coating 2610 on its distal side. This allows attachment of the protective cap 2600 to the patient's scalp.

FIG. 27 is a perspective view of a protective collar 2700 that can be disposed about a fiducial marker 2702 that has been affixed to a subject's skull. In this example, the protective collar 2700 includes a disk-like base 2704 and a circumferential cylindrical sidewall 2706 rising upward from a perimeter of the base 2704. The collar 2700 includes a radial slot 2708 in the base 2704. A first end of the radial slot 2708 terminates at an orifice 2710 at the center of the base 2704. A second end of the radial slot 2708 terminates at a peripheral slot 2712, at substantially a right angle thereto, extending up the sidewall 2706 of the collar 2700. The collar 2700 is somewhat flexible (e.g., made of plastic), and the peripheral slot 2712 and the radial slot 2708 are sized and shaped to pass the shaft 2714 of the fiducial marker 2702 through to the center orifice 2710, where it is seated. When the shaft 2714 is seated within the center orifice 2710, a height 2716 of the sidewall 2706 of the collar 2700 is greater than a height 2718 between a top of the fiducial marker 2702 and the patient's scalp 2720. When the collar 2700 has been disposed about the fiducial marker 2702, it protects the fiducial marker 2702 against a mechanical impact.

FIG. 28 is a perspective view of the collar 2700 further including a disk-like cap 2800 that fits snugly over and around the top of the collar 2700 to house and substantially enclose the fiducial marker 2702 disposed within the collar

2700. The cap **2800** is not required, but it provides additional structural strength and helps keep clean the incision through which the fiducial marker **2702** was introduced.

FIG. **29** is a perspective view of an alternate example of a collar **2900**, similar to that illustrated in FIGS. **27 – 28**, but that omits the radial slot **2708** and the peripheral slot **2712**. In this example, the protective collar **2900** includes a disk-like base **2904** and a circumferential cylindrical sidewall **2906** rising upward from a perimeter of the base **2904**. The collar **2900** includes an orifice **2910** at the center of the base **2904**. The flexible base **2904** includes small incisions **2920** extending radially from the orifice **2910** to permit the head **2922** portion of the fiducial marker **2702** (which is larger than the orifice **2910**) to pass through the orifice **2910**. The orifice **2910** is sized to accommodate the shaft **2714** portion of the fiducial marker **2702** snugly therein. When the collar **2900** is seated against the scalp **2720**, a height **2916** of the sidewall **2906** of the collar **2900** is greater than a height **2918** between a top of the fiducial marker **2702** and the scalp **2720** of the patient. When the collar **2900** has been disposed about the fiducial marker **2702**, it protects the fiducial marker **2702** against a mechanical impact, etc. The collar **2900** can also be used in conjunction with the cap **2800** illustrated in FIG. **28**, as discussed above.

FIG. **30** is a perspective view illustrating an example of a headband **3000** for protecting fiducial markers from mechanical impact. The headband **3000** is sized and shaped to fit around the skull of a subject **3002**. The headband includes one or more fixation straps **3003**, e.g., using Velcro to attach opposing sides of the headband **3000**. In one example, the headband **3000** includes one or more pre-formed holes **3004**, which are located in relationship to each other in a manner to be suitable for placing image-guided surgical (IGS) fiducial markers at the locations **3006** of the holes when the headband **3000** is placed about the subject's head. In an alternative example, the headband **3000** does not include such holes **3004**. Instead, the user cuts holes in the headband **3000** as desired for locating the fiducial markers. In yet another example, the holes **3004** are replaced by perforation openings, so that the underlying fiducial marker only pokes through the headband as much as is needed.

FIG. **31** is a side view illustrating an example of a tubular protective guide collar **3100**. The guide collar **3100** carries a fiducial marker **3102**. The

guide collar **3100** is useful for holding and guiding the fiducial marker **3102** while it is being affixed to the patient's skull, as well as for protecting the fiducial marker **3102** after it has been affixed to the patient's skull. In this example, the tubular guide collar **3100** includes an inner diameter **3104** that is

5 large enough to receive the head **3106** of the fiducial marker **3102**. An intermediate portion of the guide collar **3100** includes a circumferential neck **3107**. The neck **3107** has a slightly smaller inner diameter than the diameter of the head **3106**. However, the neck **3107** is flexible, deformable, and/or compliant enough to pass the head **3106** through the neck **3107** when the fiducial

10 marker **3102** is affixed to the patient's skull—without pulling the fiducial marker **3102** loose from the patient's skull. This can be accomplished by constructing the guide collar **3100** of a somewhat compliant plastic, and providing appropriate neck dimensions for a particular fiducial marker head **3106**. The guide collar **3100** also optionally includes a distal flange **3108**, such as to

15 provide additional stability and to enhance vertical orientation of the guide collar **3100**. The user can hold the guide collar **3100** in place, such as by pressing two fingers against the flange **3108** to hold it against the patient's scalp. This properly holds straight and orients the fiducial marker **3102** as it is threaded into or otherwise affixed to the subject's skull. It promotes an orthogonal orientation

20 of the fiducial marker **3102** with respect to the subject's skull.

FIG. 32 is a flow chart illustrating one example of using the guide collar **3100**. At **3200**, the fiducial marker **3102** is dropped into a proximal end of the guide collar **3100**. The fiducial marker **3102** falls through the proximal tubular portion and comes to rest against the interior portion of the neck **3106**, as

25 illustrated in FIG. 31. Then, at **3202**, a distal end of the guide collar **3100** is positioned against the subject's scalp, such as by pressing down against the optional flange **3108**. At **3204**, the fiducial marker **3102** is affixed to the subject's skull, such as by inserting a screwdriver tip into the proximal end of the guide collar **3100** and into corresponding screwdriver slot(s) in the head **3106** of

30 the fiducial marker **3102**, and screwing the fiducial marker **3102** into the patient's skull. At **3206**, the guide collar **3100** can be left in place, if desired, to protect the fiducial marker **3102** against a mechanical impact. When the fiducial marker **3102** is affixed to the patient's skull, and the flange **3108** rests against the patient's scalp, the height **3110** of the guide collar **3100** is greater than the

corresponding height of the fiducial marker 3102, such that the fiducial marker head 3106 is still located within the tubular guide collar 3100. This protects the fiducial marker 3100, such as from an axial mechanical impact that otherwise might potentially drive the fiducial marker 3100 deeper into the patient's skull.

5 At 3208, the guide collar 3100 can be removed while leaving the fiducial marker 3102 affixed to the subject's skull. This can be accomplished by grasping and pulling on the guide collar 3100, or by prying under the flange 3108. As discussed above, the neck 3107 is sufficiently compliant to pass the head 3106 from the proximal portion of the hourglass-shaped guide collar 3100 to its distal
10 portion. This allows the guide collar 3100 to be removed over the top of the fiducial marker 3102 while leaving it in place. Alternatively, the fiducial marker 3102 could be affixed to the subject without using the guide collar 3100, and the guide collar 3100 could later be snapped into place over the fiducial marker 3102 to protect it against a mechanical impact, as discussed above.

15 FIG. 33 is a perspective view of an alternative guide collar 3300. In this example, the guide collar 3300 has more than one piece. In FIG. 33, the guide collar 3300 includes a cylindrical tubular guide base 3302 and a cylindrical height extender 3304. In this example, the cylindrical tubular guide base 3302 includes an optional distal flange 3306. The guide base 3302 includes a side
20 access slot 3308 that is sized and shaped to pass a shaft portion 3112 of the fiducial marker 3102. The flange 3310 includes a similar slot 3310, which is aligned with the slot 3308. The cylindrical height extender 3303 can be press-fit over the guide base 3302 snugly enough to hold these two pieces together until they are again pulled apart by the user.

25 FIG. 34 is a side sectional view of the guide base 3302, the height extender 3304, a fiducial marker 3102, and a screwdriver shaft 3400. FIG. 35 is a flow chart illustrating one example of using the guide base 3302 and the height extender 3304 of FIGS. 33 – 34. At 3500, the fiducial marker 3102 is inserted into the guide base 3302, either by dropping it in the top or by inserting its shaft
30 laterally through the side access slot 3308. At 3502, a distal portion of the guide base 3302 is placed against the subject's scalp and held in place, such as by pressing down against the optional flange 3306. At 3504, the fiducial marker 3102 is affixed to the subject's skull, such as by screwing it in such as illustrated in FIG. 34. At 3506, the height extender 3304 is slid over and snugly press-

fitted around the guide base 3302. As illustrated in FIG. 34, the height extender 3304 is taller than the affixed fiducial marker 3102. In this manner, the height extender 3304 protects the fiducial marker 3102 against a mechanical impact, such as an axial blow that might otherwise drive the fiducial marker 3102 deeper into the patient's skull. At 3508, the height extender 3304 is removed by axial pulling. At 3510, the guide base 3304 is laterally removed, thereby passing the shaft of the fiducial marker 3102 out of the slot 3308.

FIG. 36 is a schematic illustration of a subject 3600 with one or more fiducial markers 3602 affixed to the subject's skull. As discussed above, in one example, the fiducial markers 3602 include substantially spherical heads with integrated conical divot receptacles therein for mating to a remotely detectable positioning instrument. FIG. 37 illustrates schematically one example of how such fiducial marker head images 3700 appear on an image created by MR, CT, or another imaging modality. For registering the patient, it is useful to know the center locations of the fiducial marker head images 3700. However, the presence of the integrated divot may confound the fiducial marker head images 3700 somewhat.

FIG. 38 illustrates schematically one example of a template 3800 including one or more concentric rings with a center indicator (such as a bull's-eye pattern or the like) such as for assisting the user in locating the center of the fiducial marker head images 3700. In one example, the template 3800 is implemented on a physical media (e.g., a transparency) that is placed over the fiducial marker head image 3700 (e.g., on a computer display, such as the IGS workstation 400). In another example, the template 3800 is implemented by computer software (e.g., as a mouse-draggable icon or feature on a computer display, such as the IGS workstation 400) that is moved using a mouse or other computer input device to place it over a fiducial marker head image 3700. In either example, the template 3800 is concentrically aligned (e.g., using one or more of its concentric rings or similar curves for aligning with a two-dimensional image of a sphere) to one of the fiducial marker head images 3700. This provides an indication of the center of that fiducial marker head image 3700. In the physical media example, the user moves a cursor to align the fiducial marker head image 3700 with the center of the template, and clicks a mouse button to select the center of the fiducial marker head image 3700. In the

software template 3800 example, the user clicks a mouse button when the software template 3800 is aligned with a center of the fiducial marker head image 3700 to select the same. The selected center of the fiducial marker head image 3700 is then used, during the patient registration process, to correlate to
5 the physical location of the apex of the conical divot, as located by the tip of the positioning device that mates thereto, as discussed above.

Although the above examples illustrated with respect to FIGS. 11 -38 have been discussed with particular emphasis on a spherical imagable fiducial marker with integrated receptacle and bone screw, it should be understood that in
10 an alternative embodiment, such examples are implemented using a cylindrical or faceted columnar shaped fiducial marker with integrated receptacle and bone screw. Moreover, in a further example, such fiducial markers include reflective outer surfaces that are recognizable by a remote positioning system, as discussed elsewhere in this document. Still further, such fiducial markers can incorporate
15 anti-microbial properties, such as by using an anti-microbial coating, or using silver or silver-based alloys for their manufacture.

In further examples, the various above-described locators (e.g., on the subject's skull, or on a wand, as illustrated in Figure 3) alternatively or additionally include an electromagnetic (EM) coil that permits determination of
20 the position of the locator using an EM coil detecting positioning system coupled to an IGS workstation rather than the optical positioning system 412 discussed above.

It is to be understood that the above description is intended to be illustrative, and not restrictive. For example, the above-described embodiments
25 may be used in combination with each other. Many other embodiments will be apparent to those of skill in the art upon reviewing the above description. The scope of the invention should, therefore, be determined with reference to the appended claims, along with the full scope of equivalents to which such claims are entitled. In the appended claims, the terms "including" and "in which" are
30 used as the plain-English equivalents of the respective terms "comprising" and "wherein." Moreover, in the following claims, the terms "first," "second," and "third," etc. are used merely as labels, and are not intended to impose numerical requirements on their objects.